

FEATURES:

1. Webbing rated at 800 lbs.
2. FDA approved.
3. **Pedi-Spider™** secures to long spine board, limiting lateral movement.
4. **Pedi-Spider™** keeps child secure to board with 3 vertical & 5 horizontal straps.
5. If you use a short back board or other immobilization device to immobilize your pediatric patient, secure them to a long spine board with the **Pedi-Spider™**.
6. Comes in a convenient vented carry case.


CAUTIONS:

1. Inspect the **Pedi-Spider™** before each use.
2. Discard the strap if you observe: broken stitching, frayed or cut fabric or straps, cracked or deteriorated hook and loop and/or other signs of deterioration.
3. Do not overtighten to assure adequate circulation.
4. **NOT INTENDED FOR VERTICAL LIFTING.**
5. **DO NOT PLACE CHEST STRAP OVER PATIENT'S DIAPHRAGM.**

Federal law restricts this device for sale by or on the order of a physician or other licensed health care professional



Explanation of Symbols	
Symbol	Explanation
	Manufacturer
	European Authorized Representative
	Refer to for Additional Information
	Caution, Consult Accompanying Documents



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: +31 (0)70 345 8570
Fax: +31 (0)70 346 7299



Emergency Products & Research
890 West Main Street
Kent, OH 44240-2218
+1-330-673-5003



Manufactured Exclusively By



Creating Industry Standards.

Emergency Products & Research, Inc.
890 West Main Street
Kent, OH 44240-2218


Inside USA: 1-800-322-5725
Outside USA: +1-330-673-5003

www.epandr.com

©2006 Emergency Products & Research, Inc.

PEDI-SPIDER™

Pediatric Immobilization Strap (EP-76)



User Manual (EP-81)

(Spine Board shown for illustration purposes only, not included with Pedi-Spider™)

Manufactured Exclusively By



Creating Industry Standards.

PEDI-SPIDER™

Pediatric Immobilization Strap

1. Lay child flat on long spine board.
2. Place **Pedi-Spider™** on child, with middle vertical strap in middle of chest and top strap at armpit level.
3. Undo velcro fasteners and place through handle on board on both sides and then back around and fasten it to itself.

4. Repeat with all 5 straps. Once all 5 straps are fastened, tighten all straps slightly. You do not want more than 2 fingers under straps. But, do not over tighten, as this may cause circulation or respiratory problems.


Recommended Pediatric Applications

Weight - 10-100 lbs.

Length - 20-43 inches.

Age - Infant to 12 years old.

(Spine Board shown for illustration purposes only, not included with Pedi-Spider™)



REFOLDING:

1. Lay strap out as if on patient.
2. Fold cross straps back on themselves. Secure with hook and loop.
3. Position yourself at the foot end of the strap and grasp the right and left sides of the foot strap.
4. Pick up cross straps in sequence, half in your right hand and half in your left hand, beginning with the bottom strap. Align the center strap folds toward you.
5. When entire strapping system has been picked up, fold in half and secure with the small hook and loop fastener located on the V-Strap.
6. Store near backboard or scoop stretcher until next use.

CARE:

1. Keep hook and loop clean; the hook and loop closures will support considerable weight and "safe use" if kept clean and lint free.
2. Hand wash. For machine washing, put strap in a pillow case, use cool water and a mild detergent and the gentle cycle.

For Professional Use Only!

Manufactured Exclusively By

**EMERGENCY
Products + Research**

Creating Industry Standards.

For best results in the field,
practice using the Pedi-Spider™ until
you become thoroughly familiar with it.

**AIR DRY ONLY. DO NOT PUT
IN DRYER. DO NOT AUTOCLAVE.**